

A CONSUMER APPLIANCE FOR DOMESTIC USE

FEEDBACK ON THE USE OF WELLBOX® FOR PROFESSIONAL PURPOSES

LPG® SYSTEMS has received reports of the Wellbox® device being used for professional purposes, in particular for therapeutic purposes by health professionals.

LPG® SYSTEMS wishes to take a firm stand against this type of practice, which runs counter to the strict purpose for which the Wellbox® device has been developed, i.e., for domestic use only, as stated in the device's User Guide.

LPG® SYSTEMS cannot guarantee that the Wellbox® device will not be used for professional purposes, which could result in tissue damage and/or injuries to vulnerable patients/clients.

LPG® SYSTEMS will not be held responsible for the professional use of the Wellbox® device, whether in the aesthetic or medical field.

This type of risk has a major impact on the image of the LPG® Group, a benchmark in terms of scientifically proven innovation, backed by professionals.

In addition, the new European Regulation 2017/745 on medical devices, which came into force on 26 May 2021, shows an important development to strengthen the safety of medical devices in the interest of patients.

With obligations to set up clinical evaluations, data transparency for reporting incidents and CE marking, these new requirements reinforce the safe use of medical devices while promoting access to innovation to offer new solutions for patient care.

This change in regulations in favour of patients strengthens LPG® SYSTEMS' position with regard to its medical and non-medical devices.



THE WELLBOX® DEVICE

The Wellbox® device offers all the expertise of the LPG® brand, recognised in the professional circuit for over 35 years, and reserved exclusively for home use.

The compact Wellbox® device works gently and naturally on the sleeping cell mechanism (cellular stimulation) to care for your body and face at home.

- Body: localised curves, cellulite, sagging skin, heavy legs.
- Face: wrinkles, sagging features, loss of radiance, puffiness, dark circles, double chin.

Unlike the professional technologies of Cellu M6® and Mobilift M6®, designed for physiotherapists and/or beauticians, Wellbox® is equipped with different, smaller treatment heads that can be applied directly onto bare skin.

As the power and treatment chamber inside the heads are designed for consumer use only (not for professional use), tissue capture is limited.

This ensures that the Wellbox® device is perfectly safe to use at home by a non-professional user, making it unsuitable for professional and/or therapeutic use (treatment of adherent areas, for example).

Unlike the Cellu M6® professional device, which has a multitude of settings, the Wellbox® device only allows limited modification of the tissue capture mode (sequential or continuous) and the suction power, which reduces the customisation of the treatment.

Finally, unlike the Cellu M6® device, the Wellbox® device is not a medical device and as such does not have any therapeutic protocol guaranteeing that treatments are carried out correctly, nor does it benefit from any scientific study in this field.

As a result, its use in professional centres is not authorised.

CONTRAINDICATIONS REMINDER:

- Open wounds, recent and old scars or any serious and/or progressive disease
- Fever, infection, inflammation, acute dermatitis, pain or redness
- Epilepsy problems
- Blood clotting disorders
- Heart problems
- Tingling in the limbs
- If you are pregnant, do not use the device on your stomach/abdomen or lower back
- Phlebitis (blood clots), varicose veins, or if you have undergone an operation on the veins within the last 2 months
- If taking long-term medications (anti-inflammatory drugs, anticoagulants, antibiotics or antihistamines)
- Swollen areas or areas with rosacea

